



Account number:

Order date:

Patient reference:

Comments:

Order sheet

1

Lens design:

Hoyalux iD MyStyle V+	Hoyalux iD LifeStyle V+	Hoyalux iD WorkStyle V+	Corridor length : mm
Design code: <i>(HOYA Internal)</i>	Clarity Harmony	Close Screen Space	MSV+: Corridor length will be calculated by HOYA LSV+: Specify corridor 11mm or 14mm WSV+: Default corridor lengths are 20mm (Close) and 22mm (Screen & Space) unless specified

2

Material:	Coating:	Photochromic:	Tint:
EYVIA 1.74 EYNOA 1.67 EYAS 1.60 PNX 1.53 Plastic 1.50	VP VP BlueControl SFT Diamond Finish Diamond Finish BlueControl Diamond Finish UV Control	Sun IQ Grey Sun IQ Brown	

Please refer to HOYA Product Guide for product availabilities

3

Prescription:

Sphere:	Cyl:	Axis:	Add:	Prism:	Base/Direction:	Prism:	Base/Direction:	Prism:
R								Share L/R
L								

4

Measurement data:


Monocular PD:	Eye Point Height:	Frame Corneal Distance:	Wearer's Pantoscopic Angle:	Frame Face Form Angle:
PD right : mm	EP right : mm	FCD right : mm	WPA : °	FFFA : °
PD left : mm	EP left : mm	FCD left : mm		

5

Addition value measurement procedure:

Phoropter					Inset value
Trial frame:	Addition lens in front	Addition lens at rear	Total power	Measured Near Distance:	R inset : mm
				MND : cm	L inset : mm

6



Frame measurements:

A	Frame width	:	mm
B	Frame height	:	mm
DBL	Bridge size	:	mm
ED	Effective dia.	:	mm

7

Frame type:

Rimless	Nylon	Plastic	Metal	Inline
Other :				
Minimum edge thickness	:	mm		

Frame to follow

Please send frame and a copy of this order form to your local HOYA office for fitting.

For Hoyalux iD MyStyle V+ order, please also complete "Patient Data" form overleaf.

Please return this form to HOYA Customer Service:

HOYA LENS AUSTRALIA

F: 02 9310 3099

E: hapl_CS_NSW@hoya.com

NEW ZEALAND

F: 09 630 4182

E: hapl_CS_NZ@hoya.com





HOYALUXiD MYSTYLE V+ Patient Data

Comments:

Name :
 Reference :
 Date of Birth : Gender: Male Female

Previous lenses:

Single Vision Progressive Single Vision contact lenses
 Bifocal Progressive design/type: Progressive contact lenses
 Reading Manufacturer:
 None Corridor length: mm
 Other:

Previous prescription:

Sphere: Cyl: Axis: Add: Prism: Base/Direction: Prism: Base/Direction:
 R
 L

Satisfaction with previous design:

Very dissatisfied Dissatisfied Not Dissatisfied/not satisfied Satisfied Very satisfied

Lifestyle profile (fill in both A and B):

A) Where do you spend your day ? Indicate % you spend outdoors. %

B) Table Below:

Activity / Importance	Select respective row of percentage spent outdoors (as indicated in A) and importance (1 = least important and 5 = most important) per activity for that row only.						
	Driving 	Computer & digital devices 	Reading 	Sports & Recreation 	Shopping & Going out 	Professional Life 	Watching TV
Daily life spent outdoors equal to or more than 40%							
Daily life spent outdoors more than 30% but less than 40%							
Daily life spent outdoors less than 30%							

Please return this form to HOYA Customer Service:

HOYA LENS AUSTRALIA
 F: 02 9310 3099
 E: hapL_CS_NSW@hoya.com

NEW ZEALAND
 F: 09 630 4182
 E: hapL_CS_NZ@hoya.com

